



PRIMARY PACKAGING/CONSUMER BOX

FRONT PANEL

Brand Name

Type of gloves: Powder Free Nitrile Examination Gloves

Non-Sterile, Finger tip textured, Ambidextrous, Single Use, Beaded Cuff

Store in a cool dry place and avoid direct sunlight.

Size of gloves, Quantity of gloves by weight



Manufactured by:
Innolates (Thailand) Limited
E1-6, Export Processing Zone, Southern Industrial Estate, Village 4,
Chalung Sub-district, Hatyai District, Songkhla 90110, Thailand.
www.karex.com.my



Imported by: Importer name and address (if any)

Regulatory compliance marking and Authorized Representative



Advena Limited (Malta)
Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta.

and/ or



Advena Limited (UK)
Pure Offices, Plato Close,
Warwick CV34 6WE, United Kingdom

BACK PANEL

Type of gloves: Powder Free Nitrile Examination Gloves

Size of gloves, Quantity of gloves by weight

Non-Sterile, Latex Free, Ambidextrous, Single Use only, Powder Free

LEFT OR RIGHT SIDE PANELS

Brand Name

Type of gloves: Powder Free Nitrile Examination Gloves

Size of gloves, Quantity of gloves by weight

Powder Free



TOP PANEL

Brand Name

Type of gloves: Powder Free Nitrile Examination Gloves

Non-Sterile, Latex Free, Ambidextrous, Single Use only, Powder Free

Size of gloves, Quantity of gloves by weight

BOTTOM PANEL

Brand Name

Type of gloves: Powder Free Nitrile Examination Gloves

Non-Sterile, Latex Free, Ambidextrous, Single Use only, Powder Free

Size of gloves, Quantity of gloves by weight



Warning and/ or precautions: This product contain Dithiocarbamates which may cause allergic reaction in some people.

Keep out of reach of children.

Avoid direct sunlight, fluorescent light, heat and moisture.

LOT



YYYY/MM



YYYY/MM

Note: LOT and EXP can be used instead of symbols.

UDI (Barcode)

Made In Thailand






















INSTRUCTION FOR USE (IFU)

Donning Gloves:

1. Before donning gloves, always perform hand hygiene.
2. Take a pair of gloves from the glove packaging. Make sure to select the right glove size.
3. Place fingers into first glove, and pull the glove onto that hand.
4. With the non-gloved hand, place fingers into second glove, and pull the glove onto that hand.
5. If you notice any rips or tears in the glove, then replace the glove.
6. Remember that if you are wearing a gown, the gloves go on top of the gown.

Doffing/Removing Gloves:

1. The outside of the gloves are contaminated.
2. Using a gloved hand, grasp the palm area of the other gloved hand and remove first glove.
3. Hold the removed glove in the gloved hand.
4. Slide fingers of ungloved hand under the remaining glove at the wrist and peel off second glove.
5. Discard gloves.
6. Remember that gloves do not replace hand hygiene. You must perform hand hygiene when taking off gloves.

	Medical Devices		Not Made with Natural Rubber Latex		Keep away from sunlight
	Store in a cool, dry place		Recyclable packaging		Do not reuse/ For single use only
	Do not use if the Package is damaged		Temperature limitation		Non-Sterile
	Manufacturer		Authorized Representative in the European Union		Authorized Representative in the United Kingdom
	Lot/ Batch code		Date of Manufacturer		Use by Date
	CE Mark (in the correct font, minimum 5mm height)		UKCA Mark (in the correct font, minimum 5mm height)		Importer
	Unique Device Identifier				
<i>Remark: Add in definition of any other symbols used on the packaging</i>					